
**THUASNE®**


## SpryStep® Footplate

<b>fr</b>	<b>Semelle semi-rigide en matériaux composites</b> .....	<b>4</b>
<b>en</b>	<b>Semi-rigid composite materials footplate</b> .....	<b>6</b>
<b>de</b>	<b>Halbstarre Einlegesohle aus Verbundmaterial</b> .....	<b>8</b>
<b>nl</b>	<b>Halfstijve zool uit composietmaterialen</b> .....	<b>10</b>
<b>it</b>	<b>Plantare semi-rigido in materiali compositi</b> .....	<b>12</b>
<b>es</b>	<b>Plantilla semirrígida de materiales compuestos</b> .....	<b>13</b>


1

2


3


4


5

6

A

(a) ADULT VERSION						
7 - 7,5 cm / 2¼ - 3" (B)  < 2 cm / ¾" (C)  						
(b) Shoe size					(f) Foot length (cm)	(f) Foot length (inch)
(c) Size	US		UK	EUR		
	(d) Men	(e) Women				
S	6 - 8½	5 - 10	3,5 - 5,5	36 - 39	22 - 26	8½ - 10½
M	7 - 10	8 - 12	5,5 - 8	39 - 42	24 - 28	9½ - 11
L	9 - 12	11 - 15	8 - 10½	42 - 45	26 - 30	10 - 11⅝
XL	11 - 14	15 - 18	10½ - 13	45 - 48	27 - 32	10⅝ - 12½
XXL	13 - 16	-	13 - 15½	48 - 51	28 - 33	11 - 13
3XL*	15 - 18	-	-	51 - 54	32 - 34	12½ - 13½

(g) PEDIATRIC VERSION					
(b) Shoe size				(f) Foot length (cm)	(f) Foot length (inch)
(c) Size	US	UK	EUR		
4XS	8 - 10	7 - 9	24 - 27	15.5 - 18	6⅞ - 7⅞
3XS	10 - 12	9 - 11,5	27 - 30	16 - 20	6⅜ - 7⅙
XXS	12 - 2	11,5 - 1	30 - 33	18 - 22	7 - 8⅞
XS	2 - 4	1 - 3,5	33 - 36	19 - 24	7⅞ - 9½

	(a)	(b)	(c)	(d)	(e)	(f)	(g)
<b>fr</b>	Version adulte	Pointure	Taille	Homme	Femme	Longueur du pied	Version pédiatrique
<b>de</b>	Ausführung für Erwachsene	Schuhgröße	Größe	Männer	Damen	Fußlänge	Ausführung für Kinder
<b>nl</b>	Versie voor volwassenen	Schoenmaat	Maat	Heren	Dames	Lengte van de voet	Versie voor kinderen
<b>it</b>	Versione per adulto	Numero di scarpa	Taglia	Uomo	Donna	Lunghezza del piede	Versione pediatrica
<b>es</b>	Versión adulto	Número de pie	Talla	Hombre	Mujer	Longitud del pie	Versión pediátrica

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## **en** **SEMI-RIGID COMPOSITE MATERIALS FOOTPLATE**

### **Description/Destination**

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

This device is a dynamic semi-rigid footplate that limits excessive movements in the forefoot.

This device is available in adult and pediatric version.

This device must be worn in combination with a soft insole and only in a closed shoe.

Right or left models available.

### **Composition**

Carbon fibre - glass fibre - epoxy resin.

### **Properties/Mode of action**

The footplate, made of composite materials, is positioned under the foot to limit movement.

The blue translucent area on the footplate is trimmable.

### **Indications**

Idiopathic toe walking.

Metatarsophalangeal joint osteoarthritis (forefoot).

Midfoot osteoarthritis.

Grade I lesions of the plantar plate of the first metatarsal joint (Turf-toe).

Morton's neuroma\*.

\*For this indication, it is recommended to use the SpryStep® Footplate in combination with soft insole.

### **Contraindications**

Do not use the product if the diagnosis has not been confirmed.

Do not apply the product in direct contact with broken skin or an open wound without an adequate dressing.

Do not use in the event of known allergy to any of the components.

Do not use for patients weighing > 135 kg / 300 lbs.

Pediatric version: do not use for patients weighing > 60 kg / 130 lbs.

Severe loss of sensation in the lower limb.

Open ulcers or scarring of the foot.

Moderate to severe edema of the foot.

Mild, moderate to severe ankle and foot spasticity.

Triplanar instability.

Plantarflexion contracture.

Moderate to severe ankle and foot deformities.

Moderate to severe ankle instabilities.

Do not use in case of unstable fractures.

Do not use this product as a Chopart footplate.

### **Precautions**

Verify the product's integrity before every use.

Do not use the device if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart.

The initial fitting and adjustment must be done by a healthcare professional.

#### Pediatric version:

It is recommended that an adult supervises the application and use of the product by a child.

Regularly check the patient's growth (change in shoe size, foot to shoe size...).

It is recommended to renew the device (bigger size) if the child's shoe size (foot length) evolves by more than 2 sizes.

Strictly comply with your healthcare professional's prescription and recommendations for use. Check the condition of the affected limb and the state of the skin daily (with particular attention for patients with sensory deficit).

**CAUTION.** For diabetics, the device should only be used with a minimum of 3 mm ( 1/8") padding protection. Seek advice from a healthcare professional.

In the event of discomfort, significant hindrance, pain, variation in limb volume, abnormal sensations or change in colour of the extremities, remove the device and consult a healthcare professional.

For hygiene, security and performance reasons, do not re-use the product for another patient.

In the event of any modification in the product's performance, remove it and consult a healthcare professional.

Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...).

Do not wear the device in a medical imaging device.

The ability to drive a vehicle with the device must be assessed by a healthcare professional and according to local regulations.

**The systematic use of a sock and a soft insole for protection is recommended when wearing the device.**

Do not expose the product to extreme temperatures.

Do not heat, this device cannot be thermomolded.

Do not kneel with the device. Ⓞ

#### **Undesirable side-effect**

This device can cause skin reactions (redness, itching, burns, blisters...) or wounds of various degrees of severity. Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

#### **Instructions for use/Application**

*The device must be worn in shoes with the following features:*

- Stiff posterior heel cup sufficiently high to fully encapsulate the foot and the orthosis: recommended height 7 - 7.5 cm (2 3/4 - 3") Ⓞ ;
- Low heel height (< 2 cm (3/4"), 1.1 - 1.3 cm (≈ 1/2") higher than the forefoot area) Ⓞ ;
- Shoe closure system: self-fastening tabs or laces.

*Pediatric version:*

- Stiff posterior heel cup sufficiently high to fully encapsulate the foot and the orthosis;
- Low heel height;
- Shoe closure system: self-fastening tabs or laces.

Trainers or athletic shoes are the best type of shoes to use with the device.

*Preparation of the device by a healthcare professional:*

Ensure that the size and the model (right or left) are appropriate for the patient's needs.

Measure the foot length from the heel to the longest toe and/or ask the patient's shoe size.

Choose the appropriate size to fit the patient, referring to the size chart.

If the patient is between 2 sizes, choose the bigger product.

Make sure that the patient's footwear complies with the recommendations.

*Sole pattern and trimming:*

Note: only the blue translucent area can be cut.

Note: we recommend to wear personal protective equipment during device adjustment.

Make sure when you trim the footplate so that it matches the full template of the shoe, to avoid leaving an edge under the plantar surface.

If the shoe has a removable inner sole, use it to draw the cutting outline on the device's sole. ❶

If there is no removable inner sole, position the patient's foot to draw the cutting outline on the device's sole. ❷

Use scissors to cut out the device sole along the traced outline. ❸

Recommendation: if scissors are used to alter the device, we recommend that the edges be lightly sanded to remove any rough edges.

Use a buffing cone or fine grade sandpaper to create a smooth finish on trimmed area(s).

Mortons extension and/or reverse Mortons extension is possible but must be modified by a healthcare professional.

*Fitting the footplate:*

Take the inner sole out of the shoe, if removable.

Place the device into the shoe ❹ and then the removable insole removed previously. ❺

Ensure the device's footplate sits flat on the sole of the shoe.

Loosen the shoe laces and slide the foot inside the shoe. ❻

You may use a shoe horn if required.

*Check the fit:*

Ensure laces or other tightening mechanisms are firmly fastened.

Ensure comfort of foot, with no impingements, prior to use.

Ask the patient to walk.

#### **Care/Maintenance**

Product can be washed in accordance with the instructions shown on this leaflet and on the label.

Wash the rigid part with a moist cloth.

If the device is exposed to seawater or chlorinated water, make sure to rinse it in clear water and dry it.

#### **Storage**

Store at room temperature, preferably in the original packaging.

#### **Disposal**

Dispose of in accordance with local regulations.

#### **COMMERCIAL WARRANTY AGREEMENT AND WARRANTY LIMITATIONS**

Thuasne® offers a free, limited commercial warranty to the user, in the territory where the device was purchased, against defects in manufacturing and workmanship for a period of:

- one year for the rigid components.

The limited warranty is effective from the date of purchase of the product by the end-user.

The limited commercial warranty does not apply to any defects in manufacturing and workmanship in case of:

- any damage occurred by a usage outside the normal and intended use of the product as mentioned in the instructions for use,

- damages occurred as part of attempts to modify the product.

Any deterioration or improper trimming of the product during its modification or adjustment process by the healthcare professional when fitting device is expressly excluded from this warranty.

Any claim for this commercial warranty must be sent by the user to the entity where the product was purchased, which will forward this claim to the corresponding Thuasne® entity.

Any warranty claim will first be reviewed by Thuasne® to determine if the conditions of the limited warranty are fulfilled and do not fall into one of the cases of exclusion of the commercial warranty.

To benefit from the warranty, the buyer must mandatorily provide an original and dated proof of purchase of the product.

If the conditions of the limited warranty are fulfilled and the claim is made by the user or its legal representative (parents, guardian...) within the warranty delays indicated above, the buyer will get a new substitution product. It is expressly agreed that this commercial warranty is in addition to the legal warranties binding the entity which sold the product to the user, in accordance with the applicable local legislation in the country of purchase of the product.

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©Thuasne - 2049502 (2024-07)